



The Transparent Specialist

BURETTES




3229

BURETTES, FOR PINCHCOCK, WITH TIP CLASS B

CAT. REF	CAP. ML.	GRADUATION INTERVAL ML	TOLERANCE \pm ML.	APPROX. HEIGHT MM	Qty. Per Case	MRP / PC	Group
3229/06	10	0.05	0.05	430	8	360	F
3229/10	25	0.1	0.1	470	8	415	F
3229/13	50	0.1	0.1	720	8	430	F

3232

BURETTES, ROTAFLO, FITTED WITH ROTAFLO GP (GENERAL PURPOSE) SCREW THREAD STOPCOCKS, CLASS A



CAT. REF.	CAP. ML.	GRADUATION INTERVAL ML	TOLERANCE \pm ML	PACKING I/C/M	MRP / PC	Group
3232/06	10	0.05	0.02	1/4/-	1360	C
3232/10	25	0.1	0.05	1/4/-	1405	C
3232/13	50	0.1	0.05	1/4/-	1656	C
3232/15	100	0.2	0.1	1/4/-	2245	E

3233

BURETTES, ROTAFLO, FITTED WITH ROTAFLO GP (GENERAL PURPOSE) SCREW THREAD STOPCOCKS, CLASS B



CAT. REF.	CAP. ML.	GRADUATION INTERVAL ML	TOLERANCE \pm ML	PACKING I/C/M	MRP / PC	Group
3233/06	10	0.05	0.05	1/4/-	485	A
3233/10	25	0.1	0.1	1/4/-	450	A
3233/13	50	0.1	0.1	1/4/-	490	A
3233/15	100	0.2	0.2	1/4/-	715	E