




The Transparent Specialist

4133

CYLINDERS, GRADUATED, SINGLE METRIC SCALE, WITH POUR OUT , WITH HEXAGONAL BASE CLASS B

| CAT. REF. | CAP. ML. | GRADU- -ATION INTERVAL ML | TOLER- -ANCE ±ML | APPROX. HEIGHT MM | PACKING I/C/M | MRP / PC | Group |
|-----------|----------|------------------------------------|------------------------|-------------------------|------------------|-------------|-------|
| 4133/05 | 5 | 0.1 | 0.1 | 110 | 1/8/- | 200 | C |
| 4133/06 | 10 | 0.2 | 0.2 | 140 | 1/8/- | 225 | A |
| 4133/10 | 25 | 0.5 | 0.5 | 170 | 1/12/- | 270 | A |
| 4133/13 | 50 | 1.0 | 1.0 | 200 | 1/12/- | 295 | A |
| 4133/15 | 100 | 1.0 | 1.0 | 260 | 1/12/72 | 335 | A |
| 4133/21 | 250 | 2.0 | 2.0 | 335 | 1/8/32 | 625 | A |
| 4133/24 | 500 | 5.0 | 5.0 | 390 | 1/8/32 | 825 | A |
| 4133/28 | 1000 | 10.0 | 10.0 | 470 | 1/0/25 | 1250 | C |
| 4133/29 | 2000 | 20.0 | 20.0 | 570 | 1/0/9 | 2400 | C |




DESICCATORS

4192

DESICCATORS WITH COVER, KNOB TOP

| CAT. REF. | SIZE MM | GROUND FLANGE APPROX. I.D. MM | PACKING I/C/M | MRP / PC | Group |
|-----------|------------|--|------------------|-------------|-------|
| 4192/01 | 100mm | 105 | 0/1/- | 4500 | D |
| 4192/02 | 150mm | 154 | 1/2/12 | 4950 | D |
| 4192/03 | 200mm | 202 | 0/1/6 | 5500 | A |
| 4192/04 | 250mm | 260 | -/1/- | 8250 | A |
| 4192/05 | 300mm | 300 | -/1/- | 13950 | D |
| 4192/06 | 500mm | 500 | -/1/- | 45000 | E |



NEW




The Transparent Specialist

4220*

DESICCATORS VACUUM, WITH TABULATED COVER, STOPPER, WITH PTFE SPINDLE

| CAT. REF. | SIZE. MM | GROUND FLANGE APPROX. ID MM | INTER-CHANGABLE JOINT SIZE | PACKING I/C/M | MRP / PC | Group |
|-----------|----------|-----------------------------|----------------------------|---------------|----------|-------|
| 4220/01 | 100mm | 105 | 29/32 | 0/1/- | 9225 | D |
| 4220/02 | 150mm | 154 | 29/32 | 1/2/12 | 13300 | D |
| 4220/03 | 200mm | 202 | 29/32 | 0/1/6 | 14025 | D |
| 4220/04 | 250mm | 260 | 29/32 | -1/- | 22500 | D |
| 4220/05 | 300mm | 300 | 29/32 | -1/- | 25625 | D |
| 4220/06 | 500mm | 500 | 34/35 | -1/- | 75000 | E |



*Vacuum type Desiccators, when ground surface is cleaned & greased, will hold a vacuum of 500 mm of Mercury (Hg) over a 24 hours period.

DISHES

4251

DISHES, CRYSTALLIZING

| CAT. REF. | SIZE MM | APPROX. O.D. MM | APPROX. HEIGHT MM | QTY PER CASE | MRP/ PC | Group |
|-----------|---------|-----------------|-------------------|--------------|---------|-------|
| 4251/1708 | 100 | 100 | 50 | 16 | 330 | D |
| 4251/2109 | 150 | 150 | 75 | 12 | 575 | D |
| 4251/2310 | 190 | 190 | 100 | 2 | 970 | D |



4271

DISHES, CULTURE, PETRI

| CAT. REF. | SIZE MM | APPROX. O.D. MM | APPROX. HEIGHT MM | PACKING I/C/M | MRP/ PC | Group |
|-----------|---------|-----------------|-------------------|---------------|---------|-------|
| 4271/0803 | 50 | 50 | 17 | 0/48/- | 130 | E |
| 4271/1303 | 80 | 80 | 17 | 0/48/- | 165 | A |
| 4271/1703 | 100 | 100 | 17 | 0/48/- | 200 | A |
| 4271/2104 | 150 | 150 | 20 | 0/20/- | 495 | E |
| 4271/2404 | 200 | 200 | 20 | 0/8/- | 1075 | E |

