

# FLASKS

## 5171

### FLASKS, BOILING, FLORENCE, FLAT BOTTOM

CAT. REF.	CAP. ML.	APPROX. O.D. MM	APPROX. HEIGHT MM	APPROX. NECK O.D. MM	PACKING I/C/M	MRP/ PC	Group
5171/15	100	64	105	24	6/48/384	95	C
5171/21	250	85	140	34	6/48/192	105	A
5171/22	300	87	160	24	6/24/96	145	C
5171/24	500	105	170	34	6/24/72	150	A
5171/28	1000	131	200	40	6/12/72	222	A
5171/29	2000	167	250	50	2/8/32	415	A
5171/30	3000	185	300	50	0/2/24	1195	C
5171/32	5000	223	340	60	1/0/9	1830	C
5171/33	6000	236	340	60	1/0/9	2900	C
5171/35	10000	285	400	60	0/2/4	4200	C
5171/40	20000	350	505	70	0/1/4	7000	C



## 5211

### FLASKS, BOILING, FLAT BOTTOM, SHORT NECK, INTERCHANGEABLE JOINT.

CAT. REF.	CAP. ML.	INTERCH-ANGEABLE JOINT	APPROX. HEIGHT MM	PACKING I/C/M	MRP/ PC	Group
5211/17	150	24/29	115	0/6/288	135	C
5211/21	250	24/29	130	0/6/192	145	A
5211/24	500	24/29	148	0/6/96	210	A






The Transparent Specialist

**5371**

**FLASKS, BOILING, ROUND BOTTOM**


CAT. REF.	CAP. ML.	APPROX. O.D. MM	APPROX. HEIGHT MM	APPROX. NECK O.D. MM	PACKING I/C/M	MRP/ PC	Group
5371/15	100	64	110	24	6/48/384	85	C
5371/21	250	85	145	34	6/48/192	95	A
5371/24	500	105	178	34	6/24/96	155	A
5371/28	1000	131	212	40	6/12/72	220	A
5371/29	2000	167	260	50	2/8/32	410	A
5371/30	3000	185	310	50	0/2/24	1155	C
5371/32	5000	223	350	60	1/0/9	1800	C
5371/33	6000	236	350	60	1/0/9	2500	C
5371/35	10000	285	420	60	0/2/4	3800	C
5371/40	20000	350	505	70	0/1/4	6750	C



**5441**

**(VESSELS) REACTION, WIDE MOUTH, FLAT FLANGE  
100mm I.D., 150 mm O.D.**

CAT. REF.	CAP. ML.	APPROX. HEIGHT MM	PACKING I/C/M	MRP/ PC	Group
5441/28	1000	170	0/2/32	1260	D
5441/29	2000	205	1/2/12	1720	A
5441/32	5000	260	1/0/9	3500	C
5441/35	10000	320	0/2/4	7500	D
5441/40	20000	400	0/1/4	13250	D



**5442**

**LIDS, FOR FLASKS REACTION CAT. NO.5441. FLAT FLANGE AND INTERCHANGEABLE JOINT**

CAT. REF.	INTERCHANGEABLE JOINT SIZE		SIDE 5°	SIDE 10°	SIDE 15°	QTY. PER CASE	MRP/ PC	Group
	CENTER	PARALLEL						
5442/1924	19/26	19/26	24/29	19/26	-	2	2985	D
5442/1914	19/26	14/23	14/23	14/23	29/32	2	2985	D
5442/2914	29/32	14/23	14/23	14/23	29/32	2	2985	D
5442/1934	19/26	19/26	19/26	19/26	34/35	2	2985	D
5442/2434	24/29	19/26	19/26	19/26	34/35	2	2985	D

