



The Transparent Specialist

**8180
PIPETTES, SEROLOGICAL, CLASS A**

CAT. REF.	CAP. ML.	GRADUATION INTERVAL ML	TOLERANCE \pm ML	PACKING I/C/M	MRP/ PC	Group
8180/61*	0.1	0.01	0.006	1/8/-	335	E
8180/62*	0.2	0.01	0.006	1/8/-	335	E
8180/01	1.0	0.01	0.006	1/8/-	280	E
8180/01*	1.0	0.1	0.006	1/8/-	300	E
8180/03	2.0	0.02	0.01	1/8/-	285	E
8180/03*	2.0	0.1	0.01	1/8/-	310	E
8180/05	5.0	0.05	0.03	1/8/-	365	E
8180/05*	5.0	0.1	0.05	1/8/-	335	E
8180/06	10.0	0.1	0.05	1/8/-	385	E
8180/10	25.0	0.2	0.1	1/8/-	460	E

These pipettes are graduated to tip.



**8191
PIPETTES, SEROLOGICAL, CLASS B**

CAT. REF.	CAP. ML.	GRADUATION INTERVAL ML	TOLERANCE \pm ML	PACKING I/C/M	MRP/ PC	Group
8191/61*	0.1	0.01	0.01	1/8/-	110	E
8191/62*	0.2	0.01	0.01	1/8/-	105	E
8191/01*	1.0	0.1	0.01	1/8/-	95	E
8191/03	2.0	0.02	0.02	1/8/-	110	E
8191/05	5.0	0.05	0.05	1/8/-	120	E
8191/05*	5.0	0.1	0.1	1/8/-	110	E
8191/06	10.0	0.1	0.1	1/8/-	110	E
8191/10	25.0	0.2	0.2	8	190	E

These pipettes are graduated to tip.

