

8193

PIPETTES, SEROLOGICAL (WHITE MARKING)

| CAT. REF. | CAP. ML. | GRADUATION INTERVAL ML | TOLERANCE ± ML | PACKING I/C/M | MRP/ PC | Group |
|-----------|----------|------------------------|----------------|---------------|---------|-------|
| 8193/61* | 0.1 | 0.01 | 0.01 | 1/8/- | 110 | E |
| 8193/62* | 0.2 | 0.01 | 0.01 | 1/8/- | 105 | E |
| 8193/01 | 1.0 | 0.01 | 0.01 | 1/8/- | 95 | E |
| 8193/03 | 2.0 | 0.02 | 0.02 | 1/8/- | 110 | E |
| 8193/05 | 5.0 | 0.05 | 0.05 | 1/8/- | 120 | E |
| 8193/06 | 10.0 | 0.1 | 0.1 | 1/8/- | 125 | E |
| 8193/10 | 25.0 | 0.2 | 0.2 | 1/8/- | 190 | E |

These pipettes are graduated to tip.



8211

PIPETTES, TRANSFER, VOLUMETRIC, CLASS A

| CAT. REF. | CAP. ML. | TOLERANCE ± ML | PACKING I/C/M | MRP/ PC | Group |
|-----------|----------|----------------|---------------|---------|-------|
| 8211/01* | 1* | 0.007 | 1/8/- | 220 | E |
| 8211/03* | 2* | 0.01 | 1/8/- | 225 | E |
| 8211/05 | 5 | 0.015 | 1/8/- | 250 | E |
| 8211/06 | 10 | 0.02 | 1/8/- | 305 | E |
| 8211/09 | 20 | 0.03 | 1/8/- | 390 | E |
| 8211/10 | 25 | 0.03 | 1/8/- | 420 | E |
| 8211/13 | 50 | 0.05 | 1/8/- | 480 | E |
| 8211/15 | 100 | 0.08 | 1/8/- | 715 | E |

* not covered by IS specifications.



8213

PIPETTES, TRANSFER, VOLUMETRIC, CLASS B

| CAT. REF. | CAP. ML. | TOLERANCE ± ML | PACKING I/C/M | MRP/PC | Group |
|-----------|----------|----------------|---------------|--------|-------|
| 8213/01 | 1 | 0.015 | 1/8/- | 100 | E |
| 8213/03 | 2 | 0.02 | 1/8/- | 100 | E |
| 8213/05 | 5 | 0.03 | 1/8/- | 105 | E |
| 8213/06 | 10 | 0.04 | 1/8/- | 125 | E |
| 8213/09 | 20 | 0.06 | 1/8/- | 150 | E |
| 8213/10 | 25 | 0.06 | 1/8/- | 155 | E |
| 8213/13 | 50 | 0.1 | 1/8/- | 180 | E |
| 8213/15 | 100 | 0.16 | 1/8/- | 280 | E |

