


# STIRRERS

## STIRRER DRIVES

A 1400 RPM 3 Phase non-flameproof motor with built-in gear is supplied along with a flexible shaft. A motor coupling to couple the flexible shaft to motor is also provided. Other end of the flexible shaft is to be fitted into the chuck.


Cat.Ref.	HP	RPM
RSD.5*	0.5	192



## FLAMEPROOF STIRRER DRIVES

A 1400 RPM 3 Phase flameproof motor coupled with reduction gear is supplied along with a flexible shaft. Motor is mounted on a metal base frame, which is designed to install easily with a Glass Assembly.


Cat.Ref.	HP	RPM
FSD.5*	0.5	192



## SPEED REGULATORS

3 Phase non-flameproof regulators are available to control the speed of stirrer drives.

Cat.Ref.	Phase	Type
REG3	3	DIMMER
VFD3*	3	VFD



\* Flameproof VFD Type regulator is also available.