




The Transparent Specialist

6531

FLASKS, KJELDAHL, ROUND BOTTOM, LONG NECK

CAT. REF.	CAP. ML.	APPROX. O.D. MM.	APPROX. HEIGHT MM	PACKING I/C/M	MRP/ PC	Group
6531/15	100	61	210	2/12/-	155	B
6531/22	300	80	298	2/12/-	240	B
6531/24	500	97	325	2/8/48	470	A
6531/27	800	113	375	2/8/48	570	A



6541

FLASKS, KJELDAHL, ROUND BOTTOM INTERCHANGEABLE JOINT

CAT. REF.	CAP. ML.	APPROX. O.D. MM.	APPROX. HEIGHT MM	INTERCH-ANGEABLE JOINT	PACKING I/C/M	MRP/ PC	Group
6541/24	500	97	325	24/29	2/8/48	585	C
6541/27	800	113	375	24/29	2/8/48	650	C




VOLUMETRIC FLASKS

6751

VOLUMETRIC FLASKS, WITH INTERCHANGEABLE STOPPER CLASS A

CAT. REF.	CAP. ML.	TOLERANCE ±ML	INTER-CHANGEABLE STOPPER SIZE	PACKING I/C/M	MRP/ PC	Group
6751/01*	1	0.015	10/15	4	685	D
6751/03*	2	0.010	10/15	4	685	D
6751/05	5	0.025	10/15	4	640	D
6751/06	10	0.025	10/15	-4/-	555	D
6751/09	20	0.04	10/15	4	555	D
6751/10	25	0.04	10/15	-4/-	515	A
6751/13	50	0.06	10/15	2/4/128	515	A
6751/15	100	0.1	14/15	2/4/128	490	A
6751/19	200	0.15	14/15	1/8/48	655	D
6751/21	250	0.15	14/15	1/48/128	660	D
6751/24	500	0.25	19/20	1/48/128	880	D
6751/28	1000	0.4	19/20	-4/-	1275	D
6751/29	2000	0.6	24/25	-2/16	2065	D




* 1 ml and 2ml sizes are of test tube shape

6752

VOLUMETRIC FLASKS, WITH INTERCHANGEABLE STOPPER CLASS B

CAT. REF.	CAP. ML.	TOLERANCE ±ML	INTER-CHANGEABLE STOPPER	PACKING I/C/M	MRP/ PC	Group
6752/01*	1	0.02	10/15	8	235	B
6752/03*	2	0.03	10/15	8	235	B
6752/05	5	0.05	10/15	8	200	B
6752/06	10	0.05	10/15	-/8/-	185	B
6752/10	25	0.08	10/15	-/8/-	185	A
6752/13	50	0.12	10/15	4/8/256	200	A
6752/15	100	0.2	14/15	4/8/256	230	A
6752/19	200	0.3	14/15	2/16/96	265	A
6752/21	250	0.3	14/15	2/16/96	275	A
6752/24	500	0.5	19/20	2/16/96	370	A
6752/28	1000	0.8	19/20	-/4/-	530	D
6752/29	2000	1.2	24/25	-/4/16	1010	D




* 1 ml and 2 ml sizes are of test tube shape

6759

VOLUMETRIC FLASKS, WITH INTERCHANGEABLE STOPPER AMBER CLASS A

CAT. REF.	CAP. ML.	TOLERANCE ±ML	INTER-CHANGEABLE STOPPER	QTY. PER CASE	MRP/ PC	Group
6759/01*	1	0.015	10/15	4	820	F
6759/03*	2	0.010	10/15	4	840	F
6759/05	5	0.025	10/15	4	740	F
6759/06	10	0.025	10/15	-/4/-	650	F
6759/09	20	0.04	10/15	-/4/-	650	F
6759/10	25	0.04	10/15	-/4/-	585	F
6759/13	50	0.06	10/15	-/4/-	585	F
6759/15	100	0.1	14/15	-/4/-	610	F
6759/19	200	0.15	14/15	-/8/-	780	F
6759/21	250	0.15	14/15	-/8/-	790	F
6759/24	500	0.25	19/20	-/8/-	1085	F
6759/28	1000	0.4	19/20	-/4/-	1545	F
6759/29	2000	0.6	24/25	-/2/16	2435	F



* 1 ml. and 2 ml. size are of test tube shape.